


GELWHITE

TEETH WHITENING SYSTEM

MD - 666

Teeth Whitening System

Release myself , smile as much as you like



It is not enough to only have lipsticks.....

1 Year warranty on all equipment
Extended warranties available
Nationwide servicing
Training provided

Teeth Whitening System

MD-666



Seven LEDs have passed special process.



Control box



Remote control



goggle



Disposable sheath



Mouth retractor

Features:

- High efficiency goose pipe design, adjust angle at random, convenient for use.
- Combination of high power Japanese BLUE led.
- Tunable bleaching time with microprocessor-controlled timer presets.
- Digital indicator with audio feedback.
- Smart auto-select power: 100 to 240 Volts AC.50/60 Hz for global use.
- Small base size design, suitable use in all kinds dental clinic.
- High sensitive infrared remote control setting function.
- Enclosed disposable sheath, prevent cross- infection, clean and convenient, let patient feel relieved.
- The small base could adjust height.

Specification:

- Broad Spectrum: 430nm - 520nm(Host lamps)
395nm - 415nm(vice lamps)
- Light Output power: Up to 2500 mW/cm²
- 100/240 Volt, 50/60 Hz
- Weight: 4.7 Kg
- Max. Working Radius: 38 cm



Before



After