



I Light

YB 7000

The latest in IPL technology.

1 Year warranty on all equipment
Extended warranties available
Nationwide servicing
Training provided

IPL YB-7000

- Two Treatment Heads
- Sapphire Soft Key

The I-Light IPL YB-7000 uses the latest technology for the treatment of unwanted hair and skin rejuvenation. Physicians may treat patients of all skin types with an optimal blending of efficacy, safety and coverage rate - without compromise. Our compact design and user-friendly controls make it easily adaptable to any practice or clinic. The YB-7000 is manufactured to all CE safety regulations.

Treatment Principles:

The IPL principle selectively effects subcutaneous pigment or blood vessels and can simultaneously stimulate the production of collagen fibre which eventually improves skin flexibility and then reduces wrinkles. In addition, long wavelength procedures can easily destroy hair follicles to restrain the rebirth of hair.

Treatments:

- | | |
|----------------------|--------------------------|
| • Hair removal | • Skin rejuvenation |
| • Spider veins | • Pigmentation treatment |
| • Eliminate wrinkles | • Acne clearance |
| • Vascular treatment | |



PREFMED

PREFERRED MEDICAL TECHNOLOGIES

Distributed by Pref Med

Tel: +27 021 791 2870 Fax: +27 021 791 2871
www.prefmedtech.com



Department of Health License No.790/13534



I Light

YB 7000

The latest in IPL technology.

Technical Specifications

Light source	:	Intense Pulse Light (3 times atmospheres press)
Energy density	:	10~50J/Cm2
Spectrum range	:	520~1200nm
	:	Option: 420~1200nm 560~1200nm 690~1200nm 755~1200nm
Cooling temperature	:	< -4 – 0 Degrees C
Light spot area	:	8x40mm2
Pulse sequence	:	1,2,3,4,5 pulse
Pulse width	:	2~15ms
Pulse intermission	:	5~60ms
Repeat frequency	:	0.3 to 1Hz
Delivery	:	Sapphire Coupling
Cooling system	:	Water cooling, forced air cooling, semi-conductor cooling
Size	:	61*50*95
Power supply	:	220V / 50Hz or 110v / 60Hz
Power	:	1000W

Before And After I-light / IPL



PREFMED

PREFERRED MEDICAL TECHNOLOGIES

Distributed by Pref Med

Tel: +27 021 791 2870 Fax: +27 021 791 2871
www.prefmedtech.com



Department of Health License No.790/13534